



# Protected PCI with Impella® Guidance Protocol Example




# Protected PCI with Impella Guidance Protocol Example




University Hospitals  
Cleveland Medical Center

# Protected PCI with Impella Guidance Protocol Example



# Hemodynamic Effects of Impella® Devices



\* Under study

- References:**
1. Fincke R, et al. *J Am Coll Cardiol* 2004
  2. den Uijl CA, et al. *Eur Heart J* 2010
  3. Mendicino DD, et al. *Am Heart J* 2007
  4. Torgersen G, et al. *Crit Care* 2009
  5. Torre-Amione G, et al. *J Card Fail* 2009
  6. Suga H, et al. *Am J Physiol* 1979
  7. Suga H, et al. *Am J Physiol* 1981
  8. Burkhoff D, et al. *Am J Physiol Heart Circ Physiol* 2005
  9. Burkhoff D, et al. *Am J Physiol Heart Circ Physiol* 2005
  10. Sauer LDC, et al. *Artif Organs* 2007
  11. Meyers B, et al. *J Am Coll Cardiol* 2003
  12. Remmelink M, et al. *Catheter Cardiovasc Interv* 2007
  13. Agel RA, et al. *J Nucl Cardiol* 2009
  14. Lom K, et al. *Clin Res Cardiol* 2009
  15. Reesink KJD, et al. *Chest* 2004
  16. Esposto M, et al. *J Am Coll Cardiol* 2018
  17. Remmelink M, et al. *Catheter Cardiovasc Interv* 2010
  18. Naidu SS. *Circulation* 2011
  19. Weber DM, et al. *Cardiac Interventions Today Supplement/Aug/Sep 2009*

## INDICATIONS

**High-Risk PCI**  
The Impella 2.5®, Impella CP® and Impella CP® with SmartAssist® Systems are temporary (≤ 6 hours) ventricular support devices indicated for use during high-risk percutaneous coronary interventions (PCI) performed in elective or urgent, hemodynamically stable patients with severe coronary artery disease, when a heart team, including a cardiac surgeon, has determined high-risk PCI is the appropriate therapeutic option. Use of the Impella 2.5, Impella CP and Impella CP with SmartAssist Systems in these patients may prevent hemodynamic instability, which can result from repeat episodes of reversible myocardial ischemia that occur during planned temporary coronary occlusions and may reduce peri- and post-procedural adverse events.

**Cardiogenic Shock**  
The Impella 2.5®, Impella CP®, Impella CP® with SmartAssist®, Impella 5.0®, Impella 5.5® with SmartAssist® and Impella LP® Catheters, in conjunction with the Automated Impella Controller™ (collectively, “Impella® System Therapy”), are temporary ventricular support devices intended for short term use (≤ 4 days for the Impella 2.5, Impella CP, and the Impella CP with SmartAssist, and ≤ 14 days for the Impella 5.0, Impella 5.5 with SmartAssist and Impella LP) and indicated for the treatment of ongoing cardiogenic shock that occurs immediately (< 48 hours) following acute myocardial infarction or open heart surgery or in the setting of cardiomyopathy, including peripartum cardiomyopathy, or myocarditis as a result of isolated left ventricular failure that is not responsive to optimal medical management and conventional treatment measures (including volume loading and use of pressors and inotropes, with or without IABP). The intent of Impella System Therapy is to reduce ventricular work and to provide the circulatory support necessary to allow heart recovery and early assessment of residual myocardial function.

## Important Risk Information for Impella devices

### CONTRAINDICATIONS

The Impella 2.5, Impella CP with SmartAssist, Impella 5.0, Impella 5.5 with SmartAssist and Impella LP are contraindicated for use with patients experiencing any of the following conditions: Mural thrombus in the left ventricle; Presence of a mechanical aortic valve or heart restrictive device; Aortic valve stenosis/calcification (equivalent to an orifice area of 0.6 cm<sup>2</sup> or less); Moderate to severe aortic insufficiency (echocardiographic assessment graded as ≥ 2); Severe peripheral arterial disease (including placement of the Impella System; Significant right heart failure); Combined cardiorespiratory failure\*; Presence of an Atrial or Ventricular Septal Defect (including post-infarct VSD)\*; Left ventricular rupture\*; Cardiac tamponade\*

### \* This condition is a contraindication for the cardiogenic shock indication only.

### POTENTIAL ADVERSE EVENTS

Acute renal dysfunction, Aortic valve injury, Bleeding, Cardiogenic shock, Cerebral vascular accident, Stroke, Death, Hemolysis, Limb ischemia, Myocardial infarction, Renal failure, Thrombocytopenia and Vascular injury

## Right Side Support

The Impella RP® System is indicated for providing temporary right ventricular support for up to 14 days in patients with a body surface area ≥ 1.5 m<sup>2</sup>, who develop acute right heart failure or decompensation following left ventricular assist device implantation, myocardial infarction, heart transplant, or open-heart surgery.

## Important Risk Information for Impella devices

### CONTRAINDICATIONS

The Impella RP System is contraindicated for patients with the following conditions: Disorders of the pulmonary artery wall that would preclude placement or correct positioning of the Impella RP device; Mechanical valves; severe valvular stenosis or valvular regurgitation of the tricuspid or pulmonary valve; Mural thrombus of the right atrium or vena cava; Anatomic conditions precluding insertion of the pump; Presence of a vena caval filter or caval interruption device, unless there is clear access from the femoral vein to the right atrium that is large enough to accommodate a 22 Fr catheter.

### CONTRAINDICATIONS

The potential adverse effects (eg, complications) associated with the use of the Impella RP System: Arrhythmia, Atrial fibrillation, Bleeding, Cardiac tamponade, Cardiogenic shock, Death, Device malfunction, Hemolysis, Hepatic failure, Insertion site infection, Perforation, Phlegmasia cerulea dolens (a severe form of deep venous thrombosis), Pulmonary valve insufficiency, Respiratory dysfunction, Sepsis, Thrombocytopenia, Thrombotic vascular (non-central nervous system) complication, Tricuspid valve injury, Vascular injury, Venous thrombosis, Ventricular fibrillation and/or tachycardia. In addition to the risks above, there are other **WARNINGS** and **PRECAUTIONS** associated with Impella devices.


Visit [www.abiomed.com/important-safety-information](http://www.abiomed.com/important-safety-information) to learn more



**24/7 Impella Clinical Support and Technical Expertise**  
**1-800-422-8666 (US)**



## Careful Risks Assessment to Optimize Patient Safety High-Risk Factors Include Anatomy, Comorbidities and Procedural Complexity



**Long-Term Mortality Risks\***

- Incomplete Revasc:** 3x increase in 5-year mortality<sup>1</sup>
- Severe Calcification:** 72% increase in 5-year mortality<sup>6</sup>
- Valvular Heart Disease:** increase in 1-year mortality<sup>7</sup>

**In-hospital Mortality Risks\***

- Surgical Ineligibility:** 6-7x increase in mortality<sup>2,3</sup>
- Chronic Kidney Disease:** 2-3x increase in mortality<sup>5</sup>
- Impaired LV (<50% EF):** 15% increase in mortality<sup>4</sup>

Riley et al. *SCAI Position Statement, May 2020* doi: 10.1002/ccd.28994

1. Farooq Vet et al. *J Am Coll Cardiol*. 2016;61(19):282-294
2. Waijdo et al. *Circulation* 2014;130:2295-2301
3. McQuiry et al. *JACC Cardiovasc Interv*. 2011;4(9):1020-1027
4. Ye Z et al. *Biomed Res Int*. 2018; 87:53176. doi.org/10.1155/2018/8753176

5. Chonryen et al. *Am Heart J*. 2006;61:52:556-564
6. Bourantas et al. *Catheter Cardiovasc Interv*. 2015;85(2):199-206
7. Stefanini et al. *Eur Heart J*. 2014;35(8):2530-2540

\*As referenced in SCAI Statement